

Res-Q-Matt

Compact vacuum mattress



Res-Q-Matt vacuum mattresses are devices primarily intended to ensure immobilization of the whole body during patient transport.

Specific features

- Extremely resistant fabric
- Removable cover to make the cleaning procedures easier

Technical data

Length	2060 ± 20 mm
Width	800 ± 10 mm
Number of handles	8
Materials	PVC, Nylon, Polystyrene
Weight (not including pump)	7,8 ± 0,5 kg
Load capacity	150 kg

Accessories

QM22109A	ABS double effect pump
QM22120A	Aluminum Pmp
QM22170A	Carry bag for vacuum mattresses
QM22199A	Repair kit
ST00593A	K BELT 3 – Three belts STX 597

Class I MD compliant with UE Reg. 2017/745

RES-Q-MATT - COMPACT VACUUM
MATT. 8 HANDLES YELLOW
QM22208A



CND Classification	V0804
Registration number	102654

Rev.0 (27/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[.spencer.it](http://www.spencer.it)