

# Res-Q-Matt

## Compact vacuum mattress



Vacuum mattresses are devices for the immobilization of trauma patients. The vacuum system allows different levels of immobilization. The more rigid helps the right positioning of the patient; softer ones allow to mitigate the low-frequency vibrations typically found on emergency vehicles during transport.

### Specific characteristics

- Extremely resistant fabric
- Removable cover to make the cleaning procedures easier
- The number of handles allow better weight distribution

### Technical data

Length	2050
Width	800 ± 10mm
N° of handles	12
Materials	Al, PVC, Nylon
Weight	5 ± 0,5 kg
Operating temp.	From -30 to +70°C

### Optionals

QM22109A	ABS double effect pump
QM22120A	Aluminum Pump
QM22170A	Carry bag for vacuum mattresses
QM22199A	Repair kit
ST00593A	K BELT 3 – Three belts STX597

### Class I Med.Device according to Dir. 93/42/CEE

MOD.M006 RES-Q-MATT -  
VACUUM MATTRESS 12 HANDLES  
ORANGE



**QM22214A**

**CND Classification** V0804

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)