

CDK

Tapered vacuum mattress



CDK vacuum mattresses are devices primarily intended to ensure immobilization of the whole body during patient transport.

Specific features

Robust, non-absorbent and self-extinguishing fabric
 Perimetral rope with soft handles for easy transport
 Tapered shape to improve adherence to the patient
 9 carrying handles make it easier to transport and allow better adaptation in various rescue scenarios

Technical data

Length	2150 ± 30 mm
Width	900 ± 30 mm
Number of handles	9
Materials	PP, PVC, Nylon, Brass, Polystyrene
Weight (not including pump)	4,1 ± 0,3 kg
Load capacity	150 kg

Accessories

ST00593A	K Belt 3 – 3 belts 1 piece with bag
QM22109A	ABS double effect pump
QM22120A	Aluminum pump
QM22170A	Carry bag for vacuum mattresses
QM22199A	Repair kit

Class I MD compliant with UE Reg. 2017/745

CDK ORANGE - TAPERED VACUUM MATTRESS

QM22300A

CND Classification V0804

Registration number 102659



Rev.0 (27/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it