

RES-Q-SPLINT KIT

Vacuum splints – Kit 3 sizes with pump and bag



Res-Q-Splint splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

Specific features

- PVC sheet easy to clean
- Separate channels avoid movement and accumulation of filling material, ensuring its better distribution making more effective the immobilization of the limb
- Fast closure thanks to the integrated tear-off straps
- Non-return valve characterized by large air passages allow to reach the desired vacuum level faster
- Carry bag with handles and housing for pump on the rear side
- The supplied pump is characterized by small size and high efficiency

Technical data

Dimensions size S	560 x 460 ± 10 mm
Dimensions size M	490 x 380 x 700 ± 10 mm
Dimensions size L	750 x 540 x 900 ± 10 mm
Weight size S	295 ± 50 g
Weight size M	630 ± 50 g
Weight size L	1050 ± 50 g
Bag dimensions	550 x 410 x 160 ± 20 mm
Materials	PVC, Al, PS
Weight of the kit	3,45 ± 0,2 kg

QM22507A

MOD.Y001 RES-Q-SPLINT KIT 3 SIZES W/PUMP AND BAG

SET consisting of the following MD compliant with UE Reg. 2017/745

QM22506A – RES-Q-SPLINT SIZE S - VAC.SPL. WITH SEPAR. CHAMB.	Reg. n. 2062063	CND V0804
QM22505A – RES-Q-SPLINT SIZE M - VAC.SPL. WITH SEPAR. CHAMB.	Reg. n. 2062062	CND V0804
QM22504A – RES-Q-SPLINT SIZE L - VAC.SPL. WITH SEPAR. CHAMB.	Reg. n. 2062061	CND V0804

Rev.0 (04/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it

Rev.0 (04/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it