

RES-Q-SPLINT

Vacuum splint



Res-Q-Splint splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

Specific features

PVC sheet easy to clean

3 separate channels, limit displacement and accumulations of the filling material, ensuring a better distribution and making the immobilization more effective

Fast closure thanks to the integrated tear-off straps

Non-return valve characterized by large air passages allow to reach the desired vacuum level faster

Technical data

Longer side	490 ± 10 mm
Shorter side	380 ± 10 mm
Length	700 ± 10 mm
Material	PVC, Al, PS
Weight	630 ± 50 g

Class I MD compliant with UE Reg. 2017/745

RES-Q-SPLINT SIZE M - VAC.

SPL. WITH SEPAR. CHAMB.

QM22520A

CND Classification V0804

Registration number 102924



Rev.0 (04/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it