

RES-Q-SPLINT PLUS KIT

Vacuum splints – Kit 3 sizes with pump and bag



Res-Q-Splint Plus splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

Specific features

- PVC sheet easy to clean
- Fast closure thanks to the integrated tear-off straps
- Non-return valve characterized by large air passages allow to reach the desired vacuum level faster
- Carry bag with handles and housing for pump on the rear side
- The supplied pump is characterized by small size and high efficiency

Technical data

Dimensions size S	735 x 640 ± 10 mm
Dimensions size M	1100 x 695 ± 10 mm
Dimensions size L	1300 x 720 ± 10 mm
Weight size S	600 ± 10 g
Weight size M	850 ± 10 g
Weight size L	1050 ± 10 g
Bag dimensions	700 x 500 x 250 ± 20 mm
Materials	PVC, Al, PS
Weight of the kit	3,8 ± 0,4 kg

QM22590A

RES-Q-SPLINT PLUS KIT - UNIVERSAL VACUUM SPLINTS

SET consisting of the following MD compliant with UE Reg. 2017/745

QM22593A – RES-Q-SPLINT PLUS SIZE S- UNIVERSAL VACUUM SPLIN.	Reg. n. 102927	CND V0804
QM22592A – RES-Q-SPLINT PLUS SIZE M- UNIVERSAL VACUUM SPLIN.	Reg. n. 102926	CND V0804
QM22591A – RES-Q-SPLINT PLUS SIZE L- UNIVERSAL VACUUM SPLIN.	Reg. n. 102931	CND V0804

Rev.0 (04/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it