

RES-Q-SPLINT PLUS

Vacuum splint



Res-Q-Splint Plus splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

Specific features

- PVC sheet easy to clean
- Fast closure thanks to the integrated tear-off straps
- Non-return valve characterized by large air passages allow to reach the desired vacuum level faster
- Suitable for ankle immobilization

Technical data

Dimensions	735 x 640 ± 10 mm
Materials	PVC, Al, PS
Weight	600 ± 10 g

Class I MD compliant with UE Reg. 2017/745

RES-Q-SPLINT PLUS SIZE S-
UNIVERSAL VACUUM SPLIN.

QM22593A

CND Classification V0804

Registration number 102927



Rev.0 (04/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it