

EMSoft

Vacuum splints – Kit 2 sizes with pump and bag



EMSoft splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

Specific features

- Low storage space required
- One-way valve
- Vinyl layer to make cleaning procedures easier
- High moldability
- Fast application thanks to the strap closures
- Vacuum pump and carry bag included

Technical data

	Dimensions (mm)	Weight (g)
Size S	750 x 600	585
Size L	1100 x 800	1130
Pump	Ø 36 x 310	500

Note: Dimensions tolerance ± 10mm. Weight tolerance ± 40 g.

QM22700A

EMSOFT KIT TWO SIZES - VACUUM SPLINTS

SET consisting of the following MD compliant with UE Reg. 2017/745

QM22702A – EMSOFT SIZE S - VACUUM SPLINTS	Reg. n. 102989	CND V0804
QM22701A – EMSOFT SIZE L - VACUUM SPLINTS	Reg. n. 102963	CND V0804

Rev.0 (04/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it