EMSoft Vacuum splints – Kit 3 sizes with pump and bag



EMSoft splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

- Specific features
- Low storage space required
- One-way valve
- Vinyl layer to make cleaning procedures easier
- High moldability
- Fast application thanks to the strap closures
- Vacuum pump and carry bag included

Technical data				
	Dimensions (mm)	Weight (g)		
Size S	750 x 600	585		
Size M	900 x 700	830		
Size L	1100 x 800	1130		
Pump	Ø 36 x 310	500		

Note: Dimensions tolerance \pm 10mm. Weight tolerance \pm 40 g.

QM22707A			
EMSOFT KIT THREE SIZES - VACUUM SPLINTS			
SET consisting of the following MD compliant with UE Reg. 2017/745			
QM22702A – EMSOFT SIZE S - VACUUM SPLINTS	Reg. n. 102989	CND V0804	
QM22705A – EMSOFT SIZE M - VACUUM SPLINTS	Reg. n. 102965	CND V0804	
QM22701A - EMSOFT SIZE L - VACUUM SPLINTS	Reg. n. 102963	CND V0804	

Rev.0 (04/05/2021) UNCHECKED COPY – further revisions will be available on <u>http://support.spencer.it</u> Spencer Italia S.r.I. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it