

EMSoft Pro

Vacuum splints



Immobilization splint for traumatized limbs. Used to limit the damage and pain to the patient related with rescue operations.

Specific Characteristics

- Low storage space required
- One-way valve
- High moldability
- Fast application thanks to quick release buckles
- Vacuum pump and carry bag included

Technical data

	Dimensions (mm)	Weight (g)
Size S	750x600	860
Size M	900x700	1200
Size L	1100x800	1510
Pump	Ø36 x 310	500

Note: Dimensions tolerance ± 10mm, Weight tolerance ± 40g

Class I Med.Device according to Dir. 93/42/CEE

EMSOFT PRO KIT 3SIZES YELLOW -
VACUUM SPLINTS

QM22761

CND Classification



V0804

Rev.1 (21/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it