

# COMBIMATT

## Vacuum mattress



Combimatt vacuum mattresses are devices primarily intended to ensure immobilization of the whole body during patient transport.

### Specific features

- Robust, non-absorbent and self-extinguishing fabric
- Reinforced handles to ensure a safe grab
- The integrated foam mattress increases patient comfort when Combimatt is used as a mattress for stretcher
- Possibility to insert a spine board

### Standard equipment

- STX 597 – One piece plastic belt (3 pcs)

### Accessories

- |          |                               |
|----------|-------------------------------|
| QM22199A | Repair kit                    |
| QM22120A | Aluminium pump                |
| QM22803A | Combimatt – PVC transport bag |

### Technical data

Length	2030 ± 10 mm
Width	990 ± 10 mm
Thickness (valve zone)	100 ± 10 mm
Number of handles	8
Materials	PV
Weight (not including pump)	6,4 ± 0,5 kg
Load capacity	150 kg

Class I MD compliant with UE Reg. 2017/745

COMBIMATT - DOUBLE USE  
VACUUM/TRANSP. MATTRESS



QM22800A

CND Classification V0804

Registration number 102666

Rev.0 (27/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)