

# COMBIMATT

## Vacuum mattress



Combimatt vacuum mattresses are devices primarily intended to ensure immobilization of the whole body during patient transport.

### Specific features

- Robust, non-absorbent and self-extinguishing fabric
- Reinforced handles to ensure a safe grab
- The integrated foam mattress increases patient comfort when Combimatt is used as a mattress for stretcher
- Possibility to insert a spine board

### Standard equipment

- STX 597 – One piece plastic belt (3 pcs)
- Aluminium pump
- Combimatt – PVC transport bag

### Accessories

- QM22199A Repair kit

### Technical data

Length	2030 ± 10 mm
Width	990 ± 10 mm
Thickness (valve zone)	100 ± 10 mm
Number of handles	8
Materials	PVC, PE, Nylon, Polystyrene, Brass
Weight (not including pump)	6,4 ± 0,5 kg
Load capacity	150 kg

Class I MD compliant with UE Reg. 2017/745

MOD.E017 COMBIMATT WITH BAG  
AND PUMP QMX 120  
QM22804A



CND Classification	V0804
Registration number	1229038

Rev.0 (27/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)