

COMBIMATT

Vacuum mattress



Vacuum mattresses are devices for the immobilization of trauma patients. The vacuum system allows different levels of immobilization. The more rigid helps the right positioning of the patient; softer ones allow to mitigate the low-frequency vibrations typically found on emergency vehicles during transport.

Specific features

- Robust, non-absorbent and self-extinguishing fabric
- Reinforced handles to ensure a safe grab
- The integrated foam mattress increases patient comfort when Combimatt is used as a mattress for stretcher
- Possibility to insert a spinal board

Standard equipment

ST00597A	STX 597 - ONE PIECE PLASTIC BELT (3pcs)
QM22109A	ABS double effect pump
QM22803A	Combimatt – PVC transport bag

Accessories

QM22199A	Repair kit
----------	------------

Technical data

Length	2030 ± 10mm
Width	990 ± 10mm
Thickness (valve zone)	100 ± 10mm
N° of handles	8
Longitudinal channels	5
Materials	PVC, PE, Nylon, Al
Weight	6,4 ± 0,5 kg
Loading capacity	150 kg

Class I Med.Device according to Dir. 93/42/CEE

COMBIMATT - MATTRESS WITH
BAG AND PUMP QMX 109



QM22805A

CND Classification V0804

Registration number 1478268

Rev.1 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it