

NEXUS

Anti-shrinkage vacuum mattress



Nexus vacuum mattresses are devices primarily intended to ensure immobilization of the whole body during patient transport.

Specific features

Robust, non-absorbent and self-extinguishing fabric
 Reinforced handles to ensure a safe grab
 Embedded inserts in the lower side of the mattress, offer more support to the patient and limit the retraction of the mattress during decompression procedures

Standard equipment

STX 597 – One piece plastic belt (3pcs)
 Vacuum pump mod. 120
 Transport bag

Accessories

QM22109A	ABS double effect pump
QM22120A	Aluminum pump
QM22170A	Carry bag for vacuum mattresses
QM22199A	Repair kit

Technical data

Length	2025 ± 10 mm
Width	855 ± 10 mm
Thickness (valve zone)	100 ± 10 mm
Number of handles	8
Longitudinal channels	5
Materials	PVC, Wood, Nylon, Brass, Polystyrene
Weight (not including pump)	5,6 ± 0,5 kg
Load capacity	150 kg

Class I MD compliant with UE Reg. 2017/745

NEXUS WITH PUMP AND TRANS.
 BAG-ANTI-SH.VAC.MATT.



QM22901A

CND Classification	V0804
Registration number	102670

Rev.0 (27/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it