

MINI MATT

Paediatric vacuum mattress



Mini Matt vacuum mattresses are devices primarily intended to ensure immobilization of the whole body during patient transport.

Specific features

Versatile, easy to use and to sanitize
Wide and sturdy carrying handles
The mattress is covered with a bag to protect the mattress from wear

Technical data

Length	1300 ± 30 mm
Width	600 ± 30 mm
Number of handles	4
Materials	PVC, Nylon
Weight (not including pump)	2,2 ± 0,3 kg
Load capacity	30 kg

Accessories

QM22109A	ABS double effect pump
QM22199A	Repair kit

Class I MD compliant with UE Reg. 2017/745

MINI MATT - PAEDIATRIC VACUUM MATTRESS

QM30161A

CND Classification



V0804

Rev.0 (27/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it