

# HIRVIN

## Carrying sheet



Hirvin injury transportation sheets are devices to be used to move and transport patients in a lying or seated position.

### Specific features

- Built with fire-resistant materials
- Easy sanitation thanks to the non-absorbent materials
- Nylon stitching to ensure high resistance
- The width of the tape of the handles ensures a secure and safe grip
- Requires minimal storage space

### Accessories

- CB03000A Carrying bag orange with handles

### Reference standards

UNI EN 1865-1

### Technical data

|                                   |                   |
|-----------------------------------|-------------------|
| Length                            | 1850 ± 20 mm      |
| Width                             | 760 ± 5 mm        |
| Wrapped dimensions <sup>(1)</sup> | 150 x 80 x 760 mm |
| # handles                         | 8                 |
| Strap width                       | 50 mm             |
| Length of handles                 | 170 mm            |
| Capacity                          | 150 kg            |
| Materials                         | Nylon, PVC        |
| Weight                            | 1.4 ± 0.2 kg      |

(1) The measurement refers to wrapped sheets. This is purely indicative, as the device can be folded or wrapped in numerous ways

### Class I MD compliant with UE Reg. 2017/745

HIRVIN TRANSFER SHEET 8  
HANDLES  
QM40200A



|                     |        |
|---------------------|--------|
| CND Classification  | V0899  |
| Registration number | 170944 |

Rev.0 (15/02/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)