


# HIRVIN

## High resistant carrying sheet



Hirvin injury transportation sheets are devices to be used to move and transport patients in a lying or seated position.

### Specific features

- Built with fire-resistant materials
- Easy sanitation thanks to the non-absorbent materials
- Nylon stitching to ensure high resistance
- The width of the tape of the handles ensures a secure and safe grip
- Requires minimal storage space

### Accessories

- CB03000A Carrying bag orange with handles

### Reference standards

UNI EN 1865-1

### Technical data

Lenght	1850 ± 10 mm
Width	760 ± 5 mm
Wrapped dimensions <sup>(1)</sup>	150 x 80 x 760 mm
# handles	10
Strap width	50 mm
Lenght of handles	170 mm
Capacity	300 kg
Materials	Nylon, PVC
Weight	1.5 kg

(1) The measurement refers to wrapped sheets. This is purely indicative, as the device can be folded or wrapped in numerous ways

Class I MD compliant with UE Reg. 2017/745

HIRVIN TEN HANDLES ORANGE

QM40300A

CND Classification



V0899

Registration number

170949

Rev.0 (15/02/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)