

HP HIRVIN PLUS

High resistance carrying sheet



HP Hirvin Plus injury transportation sheets are devices to be used to move and transport patients in a lying or seated position.

Specific features

- Built with fire-resistant materials
- Easy sanitation thanks to the non-absorbent materials
- Nylon stitching to ensure high resistance
- The width of the tape of the handles ensures a secure and safe grip
- Requires minimal storage space
- Dedicated compartment to store a pillow or the supplied carrying bag

Standard equipment

- CB03000A Carrying bag orange with handles

Reference standards

UNI EN 1865-1

Technical data

Length	1850 ± 20 mm
Width	760 ± 5 mm
Wrapped dimensions ⁽¹⁾	150 x 80 x 760 mm
# handles	8
Strap width	50 mm
Length of handles	170 mm
Capacity	150 kg
Materials	Nylon, PVC
Weight	1.5 ± 0.2 kg

(1) The measurement refers to wrapped sheets. This is purely indicative, as the device can be folded or wrapped in numerous ways

Class I MD compliant with UE Reg. 2017/745

HP HIRVIN PLUS YELLOW

QM40501A

CND Classification



V0899

Registration number

170954

Rev.0 (15/02/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it