

## **WOW**

## Carrying sheet with 10 handles and footrest



WOW injury transportation sheets are devices to be used to move and transport patients in a lying or seated position.

	Specific features
	Manufactured with fireproof materials
	Easy to clean thanks to the use of non-absorbent fabrics
	Reinforced nylon webbing to ensure high resistance and durability
	The width of the tape used for the handles, ensures an

The width of the tape used for the handles, ensures an easy grabbing increasing safety

On the patient side there are handles dedicated to the patient that make him more collaborative and increase safety

The special shape, offers high lateral containment during transport

Footrest adjustable in 2 positions to increase longitudinal containment

The supplied adjustable shoulder straps ensure a better distribution of the load on the rescuer

The longitudinal aluminium rods are removable through the dedicated holes and can be folded into 3 parts. When inserted, they increase the rigidity of the sheet, allowing to distribute the load over a wider surface and allowing to limit the flexion of the device.

The device also allows the transport of the patient in sitting position

Requires low storage space

Reference standards	S
---------------------	---

UNI EN 1865-1

Technical data			
Length	2560 ± 10 mm		
Upper part width	960 ± 5 mm		
Lower part width	700 ± 5 mm		
Dimensions when stored (1)	500 x 600 x 150 mm		
Dimensions of single folded rod	74 x 450 x 30 ± 5 mm		
# Handles	10		
Strap width	50 mm		
Length of upper handles	160 mm		
Length of lower handles	260 mm		
Rod length	1210 ± 10 mm		
Number of side strap adjustments	6		
Capacity	150 kg		
Weight	4,7 ± 0,2 kg		
Materials	Nylon, PVC, Al		

<sup>(1)</sup> The measurement refers to wrapped sheets. This is purely indicative, as the device can be folded or wrapped in numerous ways

Class I MD compliant with UE Reg. 2017/745

WOW - TRANSFER SHEET 10 HANDLES WITH FOOT SUPPORT

QM40600B

**CND Classification** 

CE

V0899