


HIRVIN MAX

Bariatric carrying sheet



Hirvin max is a device intended to transfer patients from the stretcher to the bed or for transport from the rescue place to the stretcher for patients who do not have spinal injuries. This particular version is characterized by increased loading capacity.

Caratteristiche specifiche

- Built with fire-resistant materials according to EN 597-1
- Easy sanitation thanks to the non-absorbent materials
- Nylon stitching to ensure high resistance
- The width of the tape of the handles, ensure a secure and safe grip
- Require minimal storage space

Technical data

Length	2750 ± 20 mm
Width	1500 ± 20 mm
Handles	6x4
Tapes width	50 mm
Handles lenght	170 mm
Loading capacity	250 kg
Materials	Nylon, PVC
Weight	2,2 ± 0,2 kg

Reference standards

UNI EN 1865-1
EN 597-1

Class I Med. Dev. According to Reg. 2017/745/UE

HIRVIN MAX - BARITRIC
CARRYING SHEET
QM40701A



Rev.0 (12/01/2022)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it