

AIR SPLINT KIT

Inflatable splints – Kit 6 sizes with bag



Air Splint splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

Specific features

- Low storage space required
- Quick application thanks to the zip closure
- Varying the inflation pressure is possible to increase or decrease the immobilization level of the limb
- Color coded zip puller for quick size identification
- Practical bag with compartments and strap closure

Technical data

Included sizes	Arm 37 cm
	Leg 40 cm
	Arm 64 cm
	Leg 67 cm
	Arm 75 cm
	Leg 85 cm
Splint material	PVC
Folded bag dimensions	370 x 220 x 60 mm
Opened bag dimensions	1600 x 370 mm
Weight	1750 ± 50 g

Rev.0 (04/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it

QS24202A

AIR SPLINT SET 6 SIZES W/BAG 6 COMPARTMENTS

SET consisting of the following MD compliant with UE Reg. 2017/745

QS24237 - AIR SPLINT - INFLATABLE SPLINT FOR ARM 37 CM	Reg. n. 103611	CND V0804
QS24240 - AIR SPLINT - INFLATABLE SPLINT FOR LEG 40 CM	Reg. n. 103636	CND V0804
QS24264 - AIR SPLINT - INFLATABLE SPLINT FOR ARM 64 CM	Reg. n. 103628	CND V0804
QS24267 - AIR SPLINT - INFLATABLE SPLINT FOR LEG 67 CM	Reg. n. 103639	CND V0804
QS24275 - AIR SPLINT - INFLATABLE SPLINT FOR ARM 75 CM	Reg. n. 103633	CND V0804
QS24285 - AIR SPLINT - INFLATABLE SPLINT FOR LEG 85 CM	Reg. n. 103644	CND V0804

Rev.0 (04/05/2021)

 UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it