

AIR SPLINT

Inflatable splint



Air Splint splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

Specific features

- Low storage space required
- Quick application thanks to the zip closure
- Varying the inflation pressure is possible to increase or decrease the immobilization level of the limb
- Color coded zip puller for quick size identification

Technical data

Length	400 mm
Zip puller	Blue
Material	PVC
Weight	182 g

Class I MD compliant with UE Reg. 2017/745

AIR SPLINT - INFLATABLE SPLINT
FOR LEG 40 CM

QS24240A



CND Classification V0804

Registration number 103636

Rev.0 (04/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it