

AIR SPLINT

Inflatable splint



Air Splint splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

Specific features

- Low storage space required
- Quick application thanks to the zip closure
- Varying the inflation pressure is possible to increase or decrease the immobilization level of the limb

Technical data

Length	640 mm
Material	PVC
Weight	184 g

Class I MD compliant with UE Reg. 2017/745

AIR SPLINT - INFLATABLE SPLINT
FOR ARM 64 CM

QS24264A



CND Classification V0804

Registration number 103628

Rev.0 (04/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it