

## AIR SPLINT

### Inflatable splint



Air Splint splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

#### Specific features

- Low storage space required
- Quick application thanks to the zip closure
- Varying the inflation pressure is possible to increase or decrease the immobilization level of the limb

#### Technical data

|          |        |
|----------|--------|
| Length   | 750 mm |
| Material | PVC    |
| Weight   | 228 g  |

Class I MD compliant with UE Reg. 2017/745

AIR SPLINT - INFLATABLE SPLINT  
FOR ARM 75 CM

QS24275A



CND Classification V0804

Registration number 103633

Rev.0 (04/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)