

AIR SPLINT KIT

Inflatable splints – Kit 4 sizes with bag



Air Splint splints are devices used to immobilise limbs in the event of suspected injuries or fractures of the limbs, to avoid additional strain during the pre-hospital phase.

Specific features

- Low storage space required
- Quick application thanks to the zip closure
- Varying the inflation pressure is possible to increase or decrease the immobilization level of the limb
- Color coded zip puller for quick size identification

Technical data

Included sizes	Arm 37 cm
	Leg 40 cm
	Arm 64 cm
	Leg 67 cm
Splint Material	PVC
Weight	850 ± 50 g

QS24300A

AIR SPLINT SET 4 SZ.W/BAG (L.67CM-A.64CM-L.40C

SET consisting of the following MD compliant with UE Reg. 2017/745

QS24237 - AIR SPLINT - INFLATABLE SPLINT FOR ARM 37 CM	Reg. n. 103611	CND V0804
QS24240 - AIR SPLINT - INFLATABLE SPLINT FOR LEG 40 CM	Reg. n. 103636	CND V0804
QS24264 - AIR SPLINT - INFLATABLE SPLINT FOR ARM 64 CM	Reg. n. 103628	CND V0804
QS24267 - AIR SPLINT - INFLATABLE SPLINT FOR LEG 67 CM	Reg. n. 103639	CND V0804

Rev.0 (04/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it