

B-VAC

Spare canister and tubes



Kit of spare parts for B-VAC suction unit, including canister, adult and pediatric tube.

Specific features

Disposable vessel with integrated filter

Technical data

Dimensions with vessel	210x170x75 mm
Adult catheter dimensions	ø16x22 cm
Pediatric catheter dimensions	ø9x22 mm
Vessel's capacity	350 ml
Viral and bacterial filtration efficiency	99,99%
Materials	PE, PVC

B-VAC WITH ADULT AND
PAEDIATRIC TUBE

RIRV001E

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it