


TUBE WITH CONNECTOR

For Serres canister



Disposable connection element Spencer suction devices Homejet, Ambujet using Serres canisters.

Specific Characteristics

Silicone tube

Perfect coupling to ensure no vacuum loss

Technical data

Tube internal diameter	8 mm
Tube external diameter	14 mm
Tube length	20 cm
Connector input/output	11,3/13,3-11,3 mm
Weight	30g

TUBE W/CONNECTION FOR SERRES VASE

HOME/AMBUJET

RISC006C

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it