

Guedel

Oropharyngeal cannula



Guedel cannula designed for resuscitation. Once inserted into the oral cavity, prevents the tongue from falling obstructing the upper airway allowing proper airflow.

Specific Characteristics

Color coded for a fast identification of the size
Rounded terminal part to limit possible tissue damage

Technical data

Materials	PVC
Length	50 mm
Weight	4 g
Ring colour	Light blue
Latex Free	

Class I Med.Device according to Dir. 93/42/CEE

GUEDEL DISPOSABLE
CANNULA 50 mm
RM10300A



CND Classification	R010102
Registration number	138018

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it