

Guedel

Oropharyngeal cannula



Guedel cannula designed for resuscitation. Once inserted into the oral cavity, prevents the tongue from falling obstructing the upper airway allowing proper airflow.

Specific Characteristics

Color coded for a fast identification of the size
Rounded terminal part to limit possible tissue damage

Technical data

| | |
|-------------|-------|
| Materials | PVC |
| Length | 60 mm |
| Weight | 5 g |
| Ring colour | Black |
| Latex Free | |

Class I Med.Device according to Dir. 93/42/CEE

GUEDEL DISPOSABLE
CANNULA 60 mm
RM10301A



| | |
|----------------------------|---------|
| CND Classification | R010102 |
| Registration number | 138019 |

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it