

Guedel

Oropharyngeal cannula



Guedel cannula designed for resuscitation. Once inserted into the oral cavity, prevents the tongue from falling obstructing the upper airway allowing proper airflow.

Specific Characteristics

Color coded for a fast identification of the size
Rounded terminal part to limit possible tissue damage

Technical data

Materials	PVC
Length	80 mm
Weight	10 g
Ring colour	Green
Latex Free	

Class I Med.Device according to Dir. 93/42/CEE

GUEDEL DISPOSABLE
CANNULA 80 mm.
RM10303A



CND Classification	R010102
Registration number	138021

Rev.1 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it