

# Guedel

## Oropharyngeal cannulae kit



Guedel cannulae kit destined to reanimation procedures. Once a cannula has been inserted into the oral cavity, it avoids the tongue obstructs the superior respiratory tracts, allowing an appropriate air afflux at the same time.

### Specific features

- Cannulae with colour code for a more rapid measure identification
- Rounded extremity to reduce possible damages to tissues
- Sterilizable with sterilizing fluid

### Technical data

Material	PVC
Box dimensions	173x124x48 mm
Weight	195 g
Included cannulae length	40 mm
	50 mm
	60 mm
	70 mm
	80 mm
	90 mm
	100 mm
	110 mm

Class I Med.Device according to Dir. 93/42/CEE

GUEDEL-KIT CANNULE 8  
MISURE CON BOX

RM10311A

CND Classification R010102

Registration number 148337



Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)