

Guedel

10 Oropharyngeal cannulas



Guedel cannula designed for resuscitation. Once inserted into the oral cavity, prevents the tongue from falling obstructing the upper airway allowing proper airflow.

Specific features

- Color coded for a fast identification of the size
- Rounded terminal part to limit possible tissue damage
- Can be sterilized with sterilizing liquid

Technical data

Material	PVC
Length	50 mm
Weight	4 g
Ring colour	Light blue
Included cannulas	10 pz
Latex Free	

Class I Med.Device according to Dir. 93/42/CEE

GUEDEL – DISP. OROPH.
AIRWAY 50 mm box.10 Pcs.

RM10320A



CND Classification R010102

Registration number 795067

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it