


D-Guedel

Oropharyngeal cannula



Combining two materials with very different physical characteristics into one moulding, The One-Piece design eliminates the separate bite block found in all conventional Guedel airways. As well as improved safety, the One-piece also incorporates a soft tip reducing the risk of trauma to the patient.

Technical data	
ISO Size	9
Length	90 mm
Colour	Yellow

Class I Med.Device according to Reg. 2017/745/UE	
D-CANNULA GUEDEL ISO	
9.0/SIZE3	
RM10506	

Distributed by Spencer Italia s.r.l.

Rev.0 (25/08/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it