

## D-Guedel

### Oropharyngeal cannula



Combining two materials with very different physical characteristics into one moulding, The One-Piece design eliminates the separate bite block found in all conventional Guedel airways. As well as improved safety, the One-piece also incorporates a soft tip reducing the risk of trauma to the patient.

Technical data	
ISO Size	10
Length	100 mm
Colour	Red

Class I Med.Device according to Reg. 2017/745/UE  
 D-CANNULA GUEDEL ISO  
 10/SIZE4  
**RM10507**



Distributed by Spencer Italia s.r.l.

Rev.0 (25/08/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. a socio unico Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)