

DGX 02

Laryngeal mask



Disposable laryngeal mask designed to keep patient's airway open during anesthesia or unconsciousness.

Specific features

Inflatable pillow
 Disposable and sterile device
Latex Free

Technical data

Weight	14 g
Size	1
Material	PVC

Class IIa Med.Device according to Dir. 93/42/CEE

DGX 02 SIZE 1 - LARYNGEAL MASK
 RM10601C



CND Classification	01020101
Registration number	957834

Distributed by Spencer Italia s.r.l.

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it