

## DGX 02

### Laryngeal mask



Disposable laryngeal mask designed to keep patient's airway open during anesthesia or unconsciousness.

#### Specific features

Inflatable pillow  
 Disposable and sterile device  
**Latex Free**

#### Technical data

Weight	61 g
Size	5
Material	PVC
Suitable for patients	Adult 70 – 100 kg

Class IIa Med.Device according to Dir. 93/42/CEE

DGX 02 SIZE 5 - LARYNGEAL MASK  
 RM10605C



CND Classification	01020101
Registration number	310903

Distributed by Spencer Italia s.r.l.

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)