

AIR CUFFED MASK

Disposable masks



Disposable face mask designed for use in artificial respiration and cardiopulmonary resuscitation.

Specific features

- Connection standard ISO 15mm
- Transparent shell to improve patient monitoring
- Color coded fixing hook
- Disposable device with a shelf life of five years
- The kit includes 10 masks

Latex Free

Technical data

Size	2
Dimensions	65x82x68 ± 2mm
Weight ^(1 mask)	32g
Material	PVC

Class IIa Med.Device according to Dir. 93/42/CEE

AIR CUFFED MASK WITHOUT
VALVE-SIZE 2 BOX 10 PCS

RM10862A

CND Classification



R03010199

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it