

# AIR CUFFED MASK

Disposable masks



Disposable face mask designed for use in artificial respiration and cardiopulmonary resuscitation.

### Specific features

- Connection standard ISO 15mm
- Transparent shell to improve patient monitoring
- Color coded fixing hook
- Disposable device with a shelf life of five years
- The kit includes 10 masks

**Latex Free**

### Technical data

Size	4
Dimensions	98x130x72±2 mm
Weight <sup>(1 mask)</sup>	54 g
Material	PVC

Class IIa Med.Device according to Dir. 93/42/CEE

AIR CUFFED MASK WITHOUT  
VALVE-SIZE 4 BOX 10 PCS  
RM10864A



CND Classification

R03010199

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)