

# AIR CUFFED MASK

Disposable masks



Disposable face mask designed for use in artificial respiration and cardiopulmonary resuscitation.

### Specific features

- Connection standard ISO 15mm
- Transparent shell to improve patient monitoring
- Color coded fixing hook
- Disposable device with a shelf life of five years
- The kit includes 10 masks

**Latex Free**

### Technical data

Size	5
Dimensions	105x138x78±2mm
Weight <sup>(1 mask)</sup>	67 g
Material	PVC

Class IIa Med.Device according to Dir. 93/42/CEE

AIR CUFFED MASK WITHOUT VALVE-SIZE 5 BOX 10 PCS

RM10865A

CND Classification



R03010199

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)