

AIR CUFFED MASK

Kit 10 disposable masks



Disposable face masks designed for use as an adjunct to artificial respiration and cardiopulmonary resuscitation.

Specific Characteristics

Connection standard ISO 15mm
 Transparent shell to improve patient monitoring
 Inflatable pillow for maximum adherence according to the shape of patient's face
 Color coded fixing hook

Latex Free

Technical data

Dimensions	50x60x45 ± 2 mm
Material	PVC
Weight	15g
The kit contains 10 masks sz. 0	

Class IIa Med.Device according to Dir. 93/42/CEE

AIR CUFFED MASK WITH
 VALVE-SIZE 0 BOX 10 PCS

RM10870A



CND Classification R0301019

Registration number 104872

Rev.3 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it