

CUFFED MASK

Disposable masks



Disposable face mask designed for use in artificial respiration and cardiopulmonary resuscitation.

Specific features

Connection standard ISO 15mm
 Transparent shell to improve patient monitoring
 Color coded fixing hook
 Disposable device with a shelf life of five years

Latex Free

Technical data

Size	1
Dimensions	57x72x48 ± 2mm
Weight	20g
Material	PVC

Class IIa Med.Device according to Dir. 93/42/CEE

CUFFED MASK WITHOUT VALVE-
 SIZE 1

RM10890A

CND Classification



R03010199

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it