

CUFFED MASK

Disposable mask



Disposable face mask designed for use in artificial respiration and cardiopulmonary resuscitation.

Specific features

- Connection standard ISO 15mm
- Transparent shell to improve patient monitoring
- Color coded fixing hook
- Disposable device with a shelf life of five years

Latex Free

Technical data

Size	4
Dimensions	98x130x72±2 mm
Weight	54 g
Material	PVC

Class IIa Med.Device according to Dir. 93/42/CEE

CUFFED MASK WITHOUT VALVE-
SIZE 4

RM10893A

CND Classification



R03010199

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it