

CUFFED MASK

Disposable mask



Disposable face mask designed for use in artificial respiration and cardiopulmonary resuscitation.

Specific features

Connection standard ISO 15mm
Transparent shell to improve patient monitoring
Inflatable pillow for maximum adherence to patient's face
Color coded fixing hook
Disposable device with a shelf life of five years

Latex Free

Technical data

Dimensions	50x60x45 ± 2mm
Weight	15g
Material	PVC

Class IIa Med.Device according to Dir. 93/42/CEE

CUFFED MASK SIZE 0 PVC MASK
COLOUR CODED WITH VALVE
RM10900A



CND Classification

R03010199

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it