

## CUFFED MASK

Disposable mask



Disposable face mask designed for use in artificial respiration and cardiopulmonary resuscitation.

### Specific features

Connection standard ISO 15mm  
 Transparent shell to improve patient monitoring  
 Inflatable pillow for maximum adherence to patient's face  
 Color coded fixing hook  
 Disposable device with a shelf life of five years

**Latex Free**

### Technical data

Dimensions	57x72x48 ± 2mm
Weight	20g
Material	PVC

Class IIa Med.Device according to Dir. 93/42/CEE

CUFFED MASK SIZE 1 PVC MASK  
 COLOUR CODED WITH VALVE  
**RM10901A**



**CND Classification**

R03010199

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)