

CUFFED MASK

Disposable mask



Disposable face mask designed for use in artificial respiration and cardiopulmonary resuscitation.

Specific features

Connection standard ISO 15mm
 Transparent shell to improve patient monitoring
 Inflatable pillow for maximum adherence to patient's face
 Color coded fixing hook
 Disposable device with a shelf life of five years

Latex Free

Technical data

| | |
|------------|----------------|
| Dimensions | 65x82x68 ± 2mm |
| Weight | 32 g |
| Material | PVC |

Class IIa Med.Device according to Dir. 93/42/CEE

CUFFED MASK SIZE 2 PVC MASK
 COLOUR CODED WITH VALVE
 RM10902A



CND Classification

R03010199

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it