

CUFFED MASK

Disposable mask



Disposable face mask designed for use in artificial respiration and cardiopulmonary resuscitation.

Specific features

- Connection standard ISO 15mm
- Transparent shell to improve patient monitoring
- Inflatable pillow for maximum adherence to patient's face
- Color coded fixing hook
- Disposable device with a shelf life of five years

Latex Free

Technical data

Dimensions	85x110x82 ± 2mm
Weight	54 g
Material	PVC

Class IIa Med.Device according to Dir. 93/42/CEE

CUFFED MASK SIZE 3 PVC MASK
COLOUR CODED WITH VALVE
RM10903A



CND Classification

R03010199

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it