

## D-CUFFED MASK

Disposable mask



Disposable face mask designed for use in artificial respiration and cardiopulmonary resuscitation.

### Specific features

Connection standard ISO 15mm  
 Transparent shell to improve patient monitoring  
 Inflatable pillow for maximum adherence to patient's face  
 Color coded fixing hook  
 Disposable device with a shelf life of five years

**Latex Free**

### Technical data

Dimensions	85x110x82 ± 2mm
Weight	54 g
Material	PVC

Class IIa Med.Device according to Dir. 93/42/CEE

D-CUFFED MASK SIZE 3 PVC MASK  
 WITH VALVE  
 RM10923A



CND Classification

R03010199

Distributed by Spencer Italia s.r.l.

Rev.0 (27/11/20)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)