

D-CUFFED MASK

Disposable mask



Disposable face mask designed for use in artificial respiration and cardiopulmonary resuscitation.

Specific features

- Connection standard ISO 15mm
- Transparent shell to improve patient monitoring
- Inflatable pillow for maximum adherence to patient's face
- Color coded fixing hook
- Disposable device with a shelf life of five years

Latex Free

Technical data

Dimensions	85x110x82 ± 2mm
Weight	54 g
Material	PVC

Class IIa Med.Device according to Dir. 93/42/CEE

D-CUFFED MASK SIZE 4 PVC MASK
WITH VALVE
RM10924A



CND Classification

R03010199

Distributed by Spencer Italia s.r.l.

Rev.0 (27/11/20)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it