

## B-Life

### Silicon face mask



Autoclavable face mask suitable for reanimation and assisted ventilation.

#### Specific features

Transparent shell improve patient monitoring  
Soft and hypoallergenic silicon pillow  
Latex Free

#### Technical data

Materials	Silicone
Dimensions	ø 76x52 mm
Max autoclave temperature	134°C
Weight	45 g

#### Optionals

BG02001A	Laxan transparent bag for manual resuscitators
----------	--

Class IIa Med.Device according to Dir. 93/42/CEE

B-LIFE SILICON MASK SIZE 2  
RM20402B



CND Classification	R03010199
--------------------	-----------

Registration number	145631
---------------------	--------

Rev.1 (08/04/2016)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. a socio unico Collecchio (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)