

## B-Life

### Silicon face mask



Autoclavable face mask suitable for reanimation and assisted ventilation.

#### Specific features

Connection standard ISO 15/22 mm  
 Transparent shell improve patient monitoring  
 Soft and hypoallergenic silicon pillow  
 Latex Free

#### Technical data

Dimensions	135x95x70 ± 5mm
Materials	Silicone, PC
Max autoclave temperature	134°C
Weight	80g

#### Optionals

BG02001A	Laxan transparent bag for manual resuscitators
----------	--

Class IIa Med.Device according to Dir. 93/42/CEE

B-LIFE - SILICONE FACE MASK sz 5



**RM20405B**

CND Classification	R03010199
--------------------	-----------

Registration number	145635
---------------------	--------

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[spencer.it](http://spencer.it)